

Product datasheet for **TP317720L**

MYBPC3 (NM_000256) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens myosin binding protein C, cardiac (MYBPC3), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217720 representing NM_000256 Red =Cloning site Green =Tags(s)

MPEPGKKPVSAFSKKPRSVVAAGSPAVFEATERAGVKVRWQRGGSDISASNKYGLATEGTRHTLTVRE
 VGPADQGSYAVIAGSSKVKFDLKVIEAEKAEPMLAPAPAPAEATGAPGEAPAPAAELGESAPSPKGSSSA
 ALNGPTPGAPDDPIGLFVMRPQDGEVTVGGSITFSARVAGASLLKPPVVKWFKGKWDLSSKVGQHLQL
 H
 DSYDRASKVYLFELHITDAQPAFTGSYRCEVSTKDKFDCSNFNLTVEAMGTGDLDLLSAFRRTSLAGGG
 RRISDSHEDTGILDFSSLLKKRDSFRTPRDSKLEAPAEEDVWEILRQAPPSEYERIAFYGVTDLRGMLK
 RLKGMRRDEKKSTAFQKKLEPAYQVSKGHKIRLTVELADHDAEVKWLKNGQEIQMSGSKYIFESIGAKRT
 LTISQCSLADDAAYQCVVGGEKCESTELFVKEPPVLITRPLEDQLVMVGQRVEFECEVSEEGAQVKWLKDG
 VELTREETFKYRFKKDQQRHHLIINEAMLEDAGHYALCTSGGQALAEIVQEKKLEVYQSIADLMVGAKD
 QAVFKCEVSDENVRGVWLKNGKELVPDSRIKVSHIGRVHKLTIDDVTPADEADYSFVPEGFACNLSAKLH
 FMEVKIDFVPRQEPPKIHLDPCGRIPDTIVVAGNKLRLDVPISGDPAPTIVWQKAITQGNKAPARPAPD
 APEDTGDSDEWVFDKLLCETEGRVRVETTKDRSIFTVEGAKEDEGVYTVTVKNPVGEDQVNLTVKVID
 VPDAPAAPKISNVGEDSCTVQWEPPAYDGGQPILGYILERKKKSYRWMRLNFDLIQELSHEARRMIEGV
 VYEMRVYAVNAIGMSRSPASQPFMPIGPPSEPTHLAVEDVSDTTVSLKWRPPERVGAGGLDGYSVEYCP
 EGCSEWVAALQGLTEHTSILVKDLPTGARLLFRVRAHNMAGPGAPVTTEPVTVQEILQRPRLQLPRHLR
 QTIQKKVGEPVNLLIPFQKGKPRPQVTWTKEGQPLAGEEVSIRNSPTDILFIRAARRVHSGTYQVTVRIE
 NMEDKATLVLQVVDKPSPPQDLRVTDAGWGLNVALEWKPPQDVGNTELWGYTVQKADKKTMEWFTVLE
 HYR
 RTHCVVPELIIGNGYFRVFSQNMVGFSDRAATTKEPVFIPRPGITYEPPNYKALDFSEAPSFTQPLVNR
 SVIAGYTAMLCCAVRGSPKPKISWFKNGLDLGEDARFRMFSKQGVLTLEIRKPCPFDDGGIYVCRATNLQG
 EARCECRLEVRVPQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK



[View online »](#)

Predicted MW:	140.6 kDa
Concentration:	>0.1 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000247
Locus ID:	4607
UniProt ID:	Q14896
RefSeq Size:	4200
Cytogenetics:	11p11.2
RefSeq ORF:	3822
Synonyms:	CMD1MM; CMH4; cMyBP-C; FHC; LVNC10; MYBP-C
Summary:	MYBPC3 encodes the cardiac isoform of myosin-binding protein C. Myosin-binding protein C is a myosin-associated protein found in the cross-bridge-bearing zone (C region) of A bands in striated muscle. MYBPC3, the cardiac isoform, is expressed exclusively in heart muscle. Regulatory phosphorylation of the cardiac isoform in vivo by cAMP-dependent protein kinase (PKA) upon adrenergic stimulation may be linked to modulation of cardiac contraction. Mutations in MYBPC3 are one cause of familial hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)

Product images:

Coomassie blue staining of purified MYBPC3 protein (Cat# [TP317720]). The protein was produced from HEK293T cells transfected with MYBPC3 cDNA clone (Cat# [RC217720]) using MegaTran 2.0 (Cat# [TT210002]).