

Product datasheet for **TP317724L**

DEFB124 (NM_001037500) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human defensin, beta 124 (DEFB124), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217724 representing NM_001037500 Red =Cloning site Green =Tags(s) MTQLLLFLVALLVLGHVPSGRSEFKRCWKGQGACQTYCTRQETYMHLCPDASLCCLSYALKPPPVPKHEY E TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	7.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001032589
Locus ID:	245937
UniProt ID:	Q8NES8
RefSeq Size:	216
Cytogenetics:	20q11.21



[View online »](#)


RefSeq ORF: 213

Synonyms: DEFB-24

Summary: Defensins are cysteine-rich cationic polypeptides that are important in the host immunologic response to invading microorganisms. This antimicrobial protein is secreted and is a member of the beta defensin protein family. Beta defensin genes are found in several clusters throughout the genome, with this gene mapping to a cluster at 20q11.1. The encoded protein may serve to enhance innate immunity in the prostate. [provided by RefSeq, Nov 2014]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified DEFB124 protein (Cat# [TP317724]). The protein was produced from HEK293T cells transfected with DEFB124 cDNA clone (Cat# [RC217724]) using MegaTran 2.0 (Cat# [TT210002]).