

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP317754

NDUF3 (NDUFAF3) (NM_199073) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens chromosome 3 open reading frame 60 (C3orf60), transcript variant 3, 20 μg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217754 protein sequence <mark>Red</mark> =Cloning site Green=Tags(s)
	MATALALRSLYRARPSLRCPPVELPWAPRRGHRLSPADDELYQRTRISLLQREAAQAMYIDSYNSRGFMI NGNRVLGPCALLPHSVVQWNVGSHQDITEDSFSLFWLLEPRIEIVVVGTGDRTERLQSQVLQAMRQRGIA VEVQDTPNACATFNFLCHEGRVTGAALIPPPGGTSLTSLGQAAQ
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP 951047</u>
Locus ID:	25915
UniProt ID:	<u>A4FU71</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	NDUF3 (NDUFAF3) (NM_199073) Human Recombinant Protein – TP317754
RefSeq Size:	1021
Cytogenetics:	3p21.31
RefSeq ORF:	555
Synonyms:	2P1; C3orf60; E3-3; MC1DN18
Summary:	This gene encodes a mitochondrial complex I assembly protein that interacts with complex I subunits. Mutations in this gene cause mitochondrial complex I deficiency, a fatal neonatal disorder of the oxidative phosphorylation system. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2009]

Product images:



Coomassie blue staining of purified NDUFAF3 protein (Cat# TP317754). The protein was produced from HEK293T cells transfected with NDUFAF3 cDNA clone (Cat# [RC217754]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US