

Product datasheet for **TP317759L**

LIGHT (TNFSF14) (NM_003807) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tumor necrosis factor (ligand) superfamily, member 14 (TNFSF14), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217759 representing NM_003807 Red =Cloning site Green =Tags(s)
	<p>MEESVVRPSVFVWDGQTDIPFTRLGRSHRRQSCSVARVGLGLLLLMGAGLAVQGWFLQLHWRLGEMVT RLPDGPAGSWEQLIQERRSHEVNPA AHLTGANSSLTGSGGPLLWETQLGLAFLRGLSYHDGALVVTKAGY YYISKVQLGGVGCPLGLASTITHGLYKRTPRYPEELELLVSQQSPCGRATSSSRVWWDSSFLGGVHLE AGEKVVVRVLDERLVRLRDGTRSYFGAFMV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	22 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003798
Locus ID:	8740



[View online »](#)

UniProt ID: [O43557](#)

RefSeq Size: 1491

Cytogenetics: 19p13.3

RefSeq ORF: 720


Synonyms: CD258; HVEML; LIGHT; LTg

Summary: The protein encoded by this gene is a member of the tumor necrosis factor (TNF) ligand family. This protein is a ligand for TNFRSF14, which is a member of the tumor necrosis factor receptor superfamily, and which is also known as a herpesvirus entry mediator (HVEM). This protein may function as a costimulatory factor for the activation of lymphoid cells and as a deterrent to infection by herpesvirus. This protein has been shown to stimulate the proliferation of T cells, and trigger apoptosis of various tumor cells. This protein is also reported to prevent tumor necrosis factor alpha mediated apoptosis in primary hepatocyte. Two alternatively spliced transcript variant encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction

Product images:



Coomassie blue staining of purified TNFSF14 protein (Cat# [TP317759]). The protein was produced from HEK293T cells transfected with TNFSF14 cDNA clone (Cat# [RC217759]) using MegaTran 2.0 (Cat# [TT210002]).