

Product datasheet for **TP317776L**

MST3 (STK24) (NM_003576) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human serine/threonine kinase 24 (STE20 homolog, yeast) (STK24), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217776 representing NM_003576 Red =Cloning site Green =Tags(s)
	<p>MDSRAQLWGLALNKRRATLPHPGGSTNLKADPEELFTKLEKIGKGSFGEVFKGIDNRTQKVAIKIIDLE EAEDEIEDIQEITVLSQCDSPIVTKYYGSYLKDTKLWIIMEYLGGSALDLEPGPLDETQIATILREI LKGLDYHLHSEKKIHRDIKAAANVLLSEHGEVKLADFGVAGQLTDTQIKRNTFVGTFFWMAPEVIKQLAYDS KADIWSLGITAIELARGEPPHSELHPMKVLFIPKNNPPTLEGNYSKPLKEFVEACLNKEPSFRPTAKEL LKHKFILRNAKKTSYLTELIDRYKRWKAEQSHDDSSSESDAETDGQASGGSDSGDWIFTIREKDPKNLE NGALQPSDLDRNKMKDIPKRPFSQCLSTIISPLFAELKEKSQACGGNLSIEELRGAIIYLAEEACPGISD TMVAQLVQRLQRYSLSGGGTSSH</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	49.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.




[View online »](#)

RefSeq:	NP_003567
Locus ID:	8428
UniProt ID:	Q9Y6E0
RefSeq Size:	4539
Cytogenetics:	13q32.2
RefSeq ORF:	1329
Synonyms:	HEL-S-95; MST3; MST3B; STE20; STK3

Summary: This gene encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cells death. There is a pseudogene for this gene on chromosome X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2013]

Protein Families: Druggable Genome, Protein Kinase

Product images:



Coomassie blue staining of purified STK24 protein (Cat# [TP317776]). The protein was produced from HEK293T cells transfected with STK24 cDNA clone (Cat# [RC217776]) using MegaTran 2.0 (Cat# [TT210002]).