

Product datasheet for TP317795

PFAS (NM_012393) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human phosphoribosylformylglycinamide synthase (PFAS), 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC217795 representing NM_012393
Red=Cloning site Green=Tags(s)

MSPVLHFYVRPSGHEGAAPGHTRRKLQGKLPQLQGVETELCYNVNWTAEALPSAEETKMLMWLFGCPDLLL
DDVARESWLLPGSNDLLLEVGPRLNFSTPTSTNIVSVCRTGLGPVDRVETTRRYRLSFAHPPSAEVEAI
ALATLHDMTEQHFPHPIQSFSPESMPEPLNGPINILGEGRLALEKANQELGLALDSWDLDFYTKRFQEL
QRNPSTVEAFDLAQSNEHSRHWFFKQQLHVDGQKLVHSLFESIMSTQESSNPNNVLKFCDNSSAIQGKE
VRFLRPEDPTRPSRFQQQQLRHVFTAETHNFPTGVCPSFGATTGTGGRIRDVQCTGRGAHVVAGTAGY
CFGNLHIPGYNLPWEDPSFQYPGNFARPLEVAIEASNGASDYGNKFGEPVLAFARSLGLQLPDGQRREW
IKPIMFSGGIGSMEADHISKEAPEPGMEVVKVGGPVYRIGVGGGAASSVQVQGDNTSDLDFGAVQRGDPE
MEQKMNRVIRACVEAPKGNPICSLHDQAGGNGNVLKELSDPAGAIYTSRFLGDPTLNALEIWGAEQ
ESNALLLRSPNRDFLTHVSARERCPACFVGTITGDRRIVLVDRECPVRRNGQGDAPPTPLPTVDLELE
WVLGKMPRKEFFLQRKPPMLQPLALPPGLSVHQALERLRLPAVASKRYLTNKVDRSVGGGLVAQQQCVGP
LQTPLADVAVALSHEELIGAATALGEQPVKSLLDPKVAARLVAEALTNLVFALVTDLRDVKCSGNMMW
AAKLPGEAALADACEAMVAVMAALGVAVDGGKDSLSMAARVGTETVRAPGSLVISAYAVCPDITATVTP
DLKHPEGRGHLLYVALSPGQHRLGGTALAQCFSQLGEHPPDLDPENLVRAFSITQGLLKDRLLCSGHV
SDGGLVTCLEMAFAGNCGLQVDVPVPRVDVLSVLFEEPLVLEVQEPDLAQLKRYRDAGLHCLLELGH
TGEAGPHAMVRVSVNGAVVLEEPVGEALRALWEETSFLDRLQAEPKVAEEERGLRERMGPSYCLPPTFP
KASVPREPGGSPRVAILREEGSNGDREMADAFHLAGFEVWDVTMQDLCSGAIGLDTFRGVAVFVGGFSYA
DVLGSAKGWAAAFTFHPRAGAELRRFRKRPDTFSLGVCNGCQLLALLGWVGGDPNEDAAEMGPDSQPAPR
GLLLRHNLSGRYESRWASVRVGGPALMLRGMGAVLPVWSAHGEGYVAFSSPELQAQIEARGLAPLHWA
DDDGNPTEQYPLNPNGSPGGVAGICSDGRHLAVMPHPERAVRPWQWAWRPPPFDTLTTSPWLQLFINAR
NWTLEGSC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag: C-Myc/DDK
Predicted MW: 144.6 kDa



[View online »](#)

Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_036525
Locus ID:	5198
UniProt ID:	O15067 , A8K9T9 , Q6P4B4
RefSeq Size:	5338
Cytogenetics:	17p13.1
RefSeq ORF:	4014
Synonyms:	FGAMS; FGAR-AT; FGARAT; GATD8; PURL
Summary:	Purines are necessary for many cellular processes, including DNA replication, transcription, and energy metabolism. Ten enzymatic steps are required to synthesize inosine monophosphate (IMP) in the de novo pathway of purine biosynthesis. The enzyme encoded by this gene catalyzes the fourth step of IMP biosynthesis. [provided by RefSeq, Jul 2008]
Protein Pathways:	Metabolic pathways, Purine metabolism

Product images:



Coomassie blue staining of purified PFAS protein (Cat# TP317795). The protein was produced from HEK293T cells transfected with PFAS cDNA clone (Cat# [RC217795]) using MegaTran 2.0 (Cat# [TT210002]).