

Product datasheet for **TP317910M**

CACNB1 (NM_199247) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calcium channel, voltage-dependent, beta 1 subunit (CACNB1), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217910 representing NM_199247 Red =Cloning site Green =Tags(s)

MVQKTSMSRGPYPPSQEIPMEVFDPSPOGKYSKRKGRFKRSDGSTSSDTSNSFVRQGSAESYTSRPSDS
DVSLEEDREALRKEAERQALAQLEKAKTKPVAFVRTNVGYNPSPGDEVVQGVAITFEPKDFLHIKEY
NNDWWIGRLVKEGCEVGFIPSPVKLDSLRLLEQKLRQNLGSSKSGDNSSSSLGDVVTGTRRPTPPASG
NEMTNLAFELDPLELEEEEAELGEQSGSAKTSVSSVTTPPPHKGKRIPIFFKTEHVPPYDVVPSMRPIILV
GPSLKGVEVTDMMQKALFDLKHFRFDGRISITRVADISLAKRSVLNPNPSKHIIERSNTRSSLAEVQSE
IERIFELARTLQLVALDADTINHPAQLSKTSLAPIIVYIKITSPKVLQRLIKSRGKSQSKHLNVQIAASE
KLAQCPPEMFDIILDENQLEDACEHLAEYLEAYWKATHPPSSTPPNPLLNRMTATAALAASPAPVSNLQV
QVLTSLRRNLGFWGGLESSQRGSVVPQEHEHAM


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	57.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_954855
Locus ID:	782
UniProt ID:	Q02641 , Q02641-2
RefSeq Size:	1847
Cytogenetics:	17q12
RefSeq ORF:	1569
Synonyms:	CAB1; CACNLB1; CCHLB1
Summary:	The protein encoded by this gene belongs to the calcium channel beta subunit family. It plays an important role in the calcium channel by modulating G protein inhibition, increasing peak calcium current, controlling the alpha-1 subunit membrane targeting and shifting the voltage dependence of activation and inactivation. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Ion Channels: Other
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM), MAPK signaling pathway

Product images:


Coomassie blue staining of purified CACNB1 protein (Cat# [TP317910]). The protein was produced from HEK293T cells transfected with CACNB1 cDNA clone (Cat# [RC217910]) using MegaTran 2.0 (Cat# [TT210002]).