

Product datasheet for **TP318080**

FGFR1OP2 (NM_015633) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human FGFR1 oncogene partner 2 (FGFR1OP2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218080 representing NM_015633 Red =Cloning site Green =Tags(s)

MSCTIEKALADAKALVERLRDHDAAESLIEQTTALNKRVEAMKQYQEEIQELNEVARHRPRSTLVMGIQ
QENRQIRELQQENKELRTSLEEHQSALELIMSKYREQMFRLLMASKKDDPGIIMKLKEQHSKIDMVHRNK
SEFFLDASRHILEAPQHGLERRHLEANQNELQAHVDQITEMAAVMRKAIEIDEQQGCKEQERIFQLEQE
NKGLREILQITRESFLNLRKDDASESTLSALVTNSDLSLRKS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	29.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_056448
Locus ID:	26127
UniProt ID:	Q9NVK5



[View online »](#)

RefSeq Size:	2306
Cytogenetics:	12p11.23
RefSeq ORF:	759
Synonyms:	HSPC123-like; WIT3.0
Summary:	May be involved in wound healing pathway.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:

Coomassie blue staining of purified FGFR1OP2 protein (Cat# TP318080). The protein was produced from HEK293T cells transfected with FGFR1OP2 cDNA clone (Cat# [RC218080]) using MegaTran 2.0 (Cat# [TT210002]).