

Product datasheet for **TP318337**

C16orf13 (METTL26) (NM_032366) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 16 open reading frame 13 (C16orf13), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218337 representing NM_032366 Red =Cloning site Green =Tags(s)
	 MLVAAAERNKDPILHVLQRQYLDPAQRGVRVLEVASGSGQHAAHFARAFPLAEWQPSDVDQRCLDSIAAT TQAQGLTNVKAPLHLDVTWGWGHEHWGGILPQSLDLLLCINMAHVSPLRCTEGLFRAAGHLLKPRALLITYG PYAINGKISPQSNVDFDLMLRCRNPEWGLRDTALLEDLGKASGLLLERMVDMPPANNKCLIFRKN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	22.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115742
Locus ID:	84326
UniProt ID:	Q96S19



[View online »](#)

RefSeq Size:	854
Cytogenetics:	16p13.3
RefSeq ORF:	612
Synonyms:	C16orf13; JFP2

Product images:

Coomassie blue staining of purified METTL26 protein (Cat# TP318337). The protein was produced from HEK293T cells transfected with METTL26 cDNA clone (Cat# [RC218337]) using MegaTran 2.0 (Cat# [TT210002]).