

Product datasheet for **TP318379L**

C16orf86 (NM_001012984) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 16 open reading frame 86 (C16orf86), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218379 representing NM_001012984 Red =Cloning site Green =Tags(s) MASAGAERRPGVQEATVVGQGQLTEEPGSAQTSECPVAGDQFLVPAHEARGTQSEDQRPAGAASESELQ E EGPKLGEERP KPHAGALEERGPRPVVSIVRPRHGPKRKPKVKSLSLPLGLRAHLKAEALPPKLPLQEEPE DSQSESPSAKQHKKAKRSLGAPVLHAVASMVSAPLETLRERKAQRLRPLYQVNYCNPENLQAGKG DGEAEVEAEALAPVPEEGGVEQLQALLPLAGELGPGALPCPSPLVTPHALAPLGEEAGEEPGGLPSL GVSDHKAEVDKSTQVDIDKMLSVCTAPLVPLSPQYK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	33.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001013002</u>
Locus ID:	388284



[View online »](#)

UniProt ID: [Q6ZW13](#)
RefSeq Size: 1287
Cytogenetics: 16q22.1
RefSeq ORF: 951

Product images:

Coomassie blue staining of purified C16orf86 protein (Cat# [TP318379]). The protein was produced from HEK293T cells transfected with C16orf86 cDNA clone (Cat# [RC218379]) using MegaTran 2.0 (Cat# [TT210002]).