

Product datasheet for **TP318386**

PAX2 (NM_003988) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human paired box 2 (PAX2), transcript variant c, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC218386 representing NM_003988
Clone or AA Sequence:	Red =Cloning site Green =Tags(s)

MDMHCKADPFPSAMHPGHGGVNLGGVFNVRPLPDVWRQIRIVELAHQGVPCDISRQLRVSHGCVSKILG
RYYETGSIKPGVIGGSKPKVATPKVVDKIAEYKRQNPTMFAWEIRDLLAEGICDNDTVPSVSSINRIIR
TKVQQPFHPTPDGAGTGVTAPGHTIVPSTASPPVSSASNDPVGSYSINGILGIPRSNGEKRRKRDVSE
SVPNGDSQSGVDSLRLKHLRADTFTQQLEALDRVFERPSYDPVQASEHIKSEQNEYSLPALTPGLDEV
KSSLSASTNPELGSNVSGTQTPVVTGRDMASSTLPGYPPHPPTGQGSYPTSTLAGMVPEAAVGPSSSL
MSKPGRKLAEVPPCVQPTGASSPATRTATPSTRPTTRLGDSATPPY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	41.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003979
Locus ID:	5076



[View online »](#)

UniProt ID:	Q02962
RefSeq Size:	4290
Cytogenetics:	10q24.31
RefSeq ORF:	1188
Synonyms:	FSGS7; PAPRS
Summary:	<p>PAX2 encodes paired box gene 2, one of many human homologues of the <i>Drosophila melanogaster</i> gene <i>prd</i>. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional suppression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2014]</p>
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified PAX2 protein (Cat# TP318386). The protein was produced from HEK293T cells transfected with PAX2 cDNA clone (Cat# [RC218386]) using MegaTran 2.0 (Cat# [TT210002]).