

Product datasheet for TP318477

SNX27 (NM_030918) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sorting nexin family member 27 (SNX27), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218477 representing NM_030918 Red =Cloning site Green =Tags(s)
	<p>MADEDGEGIHPSAPHRNNGGGGGGGGSLHCAGNGGGGGGGPRVVRIVKSESgyGFNVRGQVSEGGQLRSI NGELYAPLQHVSavlPGGAADRAGVRKGDRILEVNHVNVEGATHKQVVDLIRAGEKELILTVLSVPPHEA DNLDPSDDSLGQSFYDYTEKQAVPISVPRYKHVEQNGEKFWVYNVYMAGRQLCSKRYREFAILHQNLKRE FANFTFPRLPgKWpFSLSEQQLDARRRGLEEYLEKVCsIRVIGESDIMQEFLSESDENYNGVSDVELRVA LPDGTTVTVRVKKNSTTDQVYQAIaAKVGMDSTTVNYFALFEVISHSFVRKLAPNEFPHKLYIQNYTSAV PGTCLTIRKWLFTTEEEILLNDNDLAVTYFFHQAVDDVKKGYIKAEESYQLQKLYEQRKMVMYLNMLRT CEGYNEIIFPHCACDSRRKGHVITAISITHFKLHACTEEGQLENQVIAFEWDEMQRWDTDEEGMAFCFEY ARGEKKPRWVKIFTpYfNYMHECFERVFCELKWRKEEY</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	59.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_112180](#)

Locus ID: 81609

UniProt ID: [Q96L92](#)

RefSeq Size: 7085


Cytogenetics: 1q21.3

RefSeq ORF: 1584

Synonyms: MRT1; MY014

Summary: This gene encodes a member of the sorting nexin family, a diverse group of cytoplasmic and membrane-associated proteins involved in endocytosis of plasma membrane receptors and protein trafficking through these compartments. All members of this protein family contain a phosphoinositide binding domain (PX domain). A highly similar protein in mouse is responsible for the specific recruitment of an isoform of serotonin 5-hydroxytryptamine 4 receptor into early endosomes, suggesting the analogous role for the human protein. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified SNX27 protein (Cat# TP318477). The protein was produced from HEK293T cells transfected with SNX27 cDNA clone (Cat# [RC218477]) using MegaTran 2.0 (Cat# [TT210002]).