

Product datasheet for TP318518

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

IL1RA (IL1RN) (NM_173842) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human interleukin 1 receptor antagonist (IL1RN), transcript variant 1,

20 µg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC218518 representing NM_173842

or AA Sequence: Red=Cloning site Green=Tags(s)

MEICRGLRSHLITLLLFLFHSETICRPSGRKSSKMQAFRIWDVNQKTFYLRNNQLVAGYLQGPNVNLEEK IDVVPIEPHALFLGIHGGKMCLSCVKSGDETRLQLEAVNITDLSENRKQDKRFAFIRSDSGPTTSFESAA

CPGWFLCTAMEADQPVSLTNMPDEGVMVTKFYFQEDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 17.1 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 776214

Locus ID: 3557

UniProt ID: P18510





RefSeq Size: 1760

Cytogenetics: 2q14.1 RefSeq ORF: 531

Synonyms: DIRA; ICIL-1RA; IL-1ra; IL-1ra3; IL-1RN; IL1F3; IL1RA; IRAP; MVCD4


Summary: The protein encoded by this gene is a member of the interleukin 1 cytokine family. This

protein inhibits the activities of interleukin 1, alpha (IL1A) and interleukin 1, beta (IL1B), and modulates a variety of interleukin 1 related immune and inflammatory responses, particularly in the acute phase of infection and inflammation. This gene and five other closely related cytokine genes form a gene cluster spanning approximately 400 kb on chromosome 2. A polymorphism of this gene is reported to be associated with increased risk of osteoporotic fractures and gastric cancer. Several alternatively spliced transcript variants encoding distinct

isoforms have been reported. [provided by RefSeq, Aug 2020]

Protein Families: Druggable Genome, Secreted Protein

Product images:



Coomassie blue staining of purified IL1RN protein (Cat# TP318518). The protein was produced from HEK293T cells transfected with IL1RN cDNA clone (Cat# [RC218518]) using MegaTran 2.0 (Cat# [TT210002]).