

Product datasheet for **TP318518M**

IL1RA (IL1RN) (NM_173842) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human interleukin 1 receptor antagonist (IL1RN), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218518 representing NM_173842 Red =Cloning site Green =Tags(s)
	MEICRGLRSHLITLLLFLFHSETICRPSGRKSSKMQAFRIWDVNQKTFYLRNNQLVAGYLQGPVNVLEEK IDVWPIEPHALFLGIHGGMCLSCVKSGDETRLQLEAVNITDLSNRKQDKRFQAFIRSDSGPTTSFESAA CPGWFLCTAMEADQPVSLTNMPDEGVMVTKFYFQEDE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_776214
Locus ID:	3557
UniProt ID:	P18510



[View online »](#)

RefSeq Size: 1760

Cytogenetics: 2q14.1


RefSeq ORF: 531

Synonyms: DIRA; ICIL-1RA; IL-1ra; IL-1ra3; IL-1RN; IL1F3; IL1RA; IRAP; MVCD4

Summary: The protein encoded by this gene is a member of the interleukin 1 cytokine family. This protein inhibits the activities of interleukin 1, alpha (IL1A) and interleukin 1, beta (IL1B), and modulates a variety of interleukin 1 related immune and inflammatory responses, particularly in the acute phase of infection and inflammation. This gene and five other closely related cytokine genes form a gene cluster spanning approximately 400 kb on chromosome 2. A polymorphism of this gene is reported to be associated with increased risk of osteoporotic fractures and gastric cancer. Several alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Aug 2020]

Protein Families: Druggable Genome, Secreted Protein

Product images:



Coomassie blue staining of purified IL1RN protein (Cat# [TP318518]). The protein was produced from HEK293T cells transfected with IL1RN cDNA clone (Cat# [RC218518]) using MegaTran 2.0 (Cat# [TT210002]).