

Product datasheet for **TP318521M**

TRIM47 (NM_033452) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tripartite motif-containing 47 (TRIM47), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218521 representing NM_033452 Red =Cloning site Green =Tags(s)
	MDGSGPFSCPICLEPLREPVTLP CGHNFCLACLGALWPHRGASGAGGPGGAARCPLCQEPFPDGLQLRK N HTLSELLQLRQSGPGSGPGPAPALAPEPSAPSALPSVPEPSAPCAPEPWPAGEEPPVRCDACPEGAALPA ALSCLSCASFPAHLGPHERSPALRGHRLVPPLRRLEESLCPRHLRPLERYCRAERVCLCEACAAQEHR GHELVPLEQERALQAEQSKVLSAVEDRMDDELGAGIAQSRRTVALIKSAAVAERERSRLFADAAAALQG FQTQVLGFIEEGEAAMLGRSQGDLRRQEEQRSRLSRARQNLSQVPEADSVSFLQELLALRLALEDGCGPG PGPPRELSFTKSSQAVRAVRDMLAVACVNQWEQLRGPGGNEDGPQKLDSEADAEPQDLESTNLLSEAP R DYFLKFAYIVDLSDTADKFLQLFGTKGVKRVLCPINYPLSPTRFTHCEQVLGEGALDRGTYWVEVEIIE GWVSMGVMAEDFSPQEPYDRGLGRNAHSCCLQWNGRSFSVWFHGLEAPLPHFPSPTVGVCLLEYADRA LA FYAVRDGKMSLLRRLKASRPRRGGIPASPIDPFQSRLDSHFAGLFTHRLKPAFFLESVDAHLQIGPLKKS CISVLKRR SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	69.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



View online »

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_258411
Locus ID:	91107
UniProt ID:	Q96LD4
RefSeq Size:	2268
Cytogenetics:	17q25.1
RefSeq ORF:	1914
Synonyms:	GOA; RNF100
Summary:	E3 ubiquitin-protein ligase that mediates the ubiquitination and proteasomal degradation of CYLD.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified TRIM47 protein (Cat# [TP318521]). The protein was produced from HEK293T cells transfected with TRIM47 cDNA clone (Cat# [RC218521]) using MegaTran 2.0 (Cat# [TT210002]).