

## Product datasheet for **TP318580L**

### BCAS4 (NM\_198799) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human breast carcinoma amplified sequence 4 (BCAS4), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218580 representing NM_198799 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MQRTGGGAPRPGRNHGLPGSLRQDPVALLMLLVDADQPEPMRSGARELALFLTPEPGAEAKEVEETIEG MLLRLEEFCSLADLIRS DTSQILEENIPVLKAKLTEMRGYAKVDRLEAFVKMVGHHVAFLEADV LQAER DHGAFPQALRRWLGSAGLPSFRNKSPAPVPVTYELPTLYRTEDYFPVDAGEAQHHPRTCPRPL  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	22.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_942094</a>
Locus ID:	55653
UniProt ID:	<a href="#">Q8TDM0</a>



[View online »](#)


RefSeq Size: 1271

Cytogenetics: 20q13.13

RefSeq ORF: 609

Synonyms: CNOL

### Product images:



Coomassie blue staining of purified BCAS4 protein (Cat# [TP318580]). The protein was produced from HEK293T cells transfected with BCAS4 cDNA clone (Cat# [RC218580]) using MegaTran 2.0 (Cat# [TT210002]).