

Product datasheet for **TP318736**

SLC34A2 (NM_006424) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human solute carrier family 34 (sodium phosphate), member 2 (SLC34A2)
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>Peptide sequence encoded by RC218736 Blue=ORF Red=Cloning site Green=Tag(s)

MAPWPELGDAQPNPKYLEGAAGQQPTAPDKSKETNKTDNTEAPVTKIELLPSYSTATLIDEPTVDDP
WNLPTLQDSGIKWSERDTKGKILCFFQIGIRLILLGLFLYFFVCSLDILSSAFQLVGGKMGAGQFFSNSS
IMSNPLLGLVIGLVTVLVQSSSTSTSIVSMVSSLLTVRAAIIIMGANIGTSITNTIVALMQVGDR
SEFRRRAFAGATVHDFFNWLSVLVLLPVEVATHYLEIITQLIVESFHFKNGEDAPDLLKVITKPFTKLIV
QLDKKVISQIAMNDEKAKNKSIVKIWCKTFTNKTQINVTVPSTANCTSPSLCWTGDIQNWTMKNVTYKE
NIAKCQHIFVNFHLPDLAVGTILLILLSLLVLCGLMIVKILGSVLKGVATVIKKTINTDFPFPAWL
TGYLAILVGAGMTFIVQSSSVFTSALTPLIGIGVITIERAYPLTLGNSIGTTTTAILAALASPGNALRS
SLQIALCHFFFNISGILLWYPIFTRLPIMAKGLGNISAKYRWFVYLIIFFLIPLTVFGLSLAGW
RVLVGVGVPVVFIIILVLCRLQLSRCPRVLPKQLQNWNLPLWMRSLKPWDAVSKFTGCFQMRCCCC
CRVCCRACLLCGCPKCCRCCKCEDLEEAQEGQDVPVKAPETFDNITISREAQGEVPASDSKTECTAL
SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Recombinant protein using RC218736 also available, [TP318736](#)


Tag:	C-Myc/DDK
Predicted MW:	75.5 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006415
Locus ID:	10568
UniProt ID:	O95436
RefSeq Size:	2280
Cytogenetics:	4p15.2
RefSeq ORF:	2079
Synonyms:	NAPI-3B; NAPI-IIb; NPTIIb; PULAM
Summary:	The protein encoded by this gene is a pH-sensitive sodium-dependent phosphate transporter. Phosphate uptake is increased at lower pH. Defects in this gene are a cause of pulmonary alveolar microlithiasis. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, May 2010]
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified SLC34A2 protein (Cat# TP318736). The protein was produced from HEK293T cells transfected with SLC34A2 cDNA clone (Cat# [RC218736]) using MegaTran 2.0 (Cat# [TT210002]).