

## Product datasheet for **TP318774M**

### Phosphoglucomutase 5 (PGM5) (NM\_021965) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphoglucomutase 5 (PGM5), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	Recombinant protein was produced with TrueORF clone, RC218774.
Tag:	C-Myc/DDK
Predicted MW:	62 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_068800</a></u>
Locus ID:	5239
UniProt ID:	<u><a href="#">Q15124</a></u>
RefSeq Size:	3341
Cytogenetics:	9q21.11
RefSeq ORF:	1701
Synonyms:	PGMRP


[View online »](#)

**Summary:** Phosphoglucomutases (EC 5.2.2.2.), such as PGM5, are phosphotransferases involved in interconversion of glucose-1-phosphate and glucose-6-phosphate. PGM activity is essential in formation of carbohydrates from glucose-6-phosphate and in formation of glucose-6-phosphate from galactose and glycogen (Edwards et al., 1995 [PubMed 8586438]).[supplied by OMIM, Mar 2008]

**Protein Families:** Druggable Genome

**Product images:**



Coomassie blue staining of purified PGM5 protein (Cat# [TP318774]). The protein was produced from HEK293T cells transfected with PGM5 cDNA clone (Cat# [RC218774]) using MegaTran 2.0 (Cat# [TT210002]).