

Product datasheet for TP318775L

NLRP8 (NM_176811) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human NLR family, pyrin domain containing 8 (NLRP8), 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC218775 representing NM_176811
Red=Cloning site **Green**=Tags(s)

MSDVNPPSDTPIPFSSSTHSSHIPPWTFSCYPGSPCENGVMLYMRNVSHEELQRFKQLLLTELSTGTMP
ITWDQVETASWAEVHLLIERFPGRAWDVTSNIFAIMNCDKMCVVRREINAILPTLEPEDLNGETQV
NLEEGESGKIRRYKSNVMEKFFPIWDITWPGNQRDFFYQGVHRHEEYLPCLLLPKRPQGRQPKTVAIQG
APGIGKTLAKKVMFEWARNKFYAHKRWCIFYHFCQEVNQTTDQSFSELIEQKWPGSQDLVSKIMSKPDQ
LLLLLDGFEELTSTLIDRLEDLSEDRQKLPGSVLLSSLLSKTMLPEATLLIMIRFTSWQTCKPLLKCP
LVTLPGFNTMEKIKYFQMYFGHTEEGDQVLSFAMENTILFSMCRVPVVCWMVCSGLKQQMERNLQSC
PNATSVFVRYISSLPTRAENFSRKIHQAQLEGLCHLAADSMWHRKWWLGKEDLEEAKLDQGTGVTAFGLM
SILRRIAGEEDHYVFTLVTFQEFFAALFYVLCFPQRLKNFVLSHVNIQRLIASPRGSKSYLSHMGLFLF
GFLNEACASAVEQSFQCKVSGFNKRKLLKVIPLLHKCDPPSPGSGVPQLFYCLHEIREEAFVSQALNDYH
KVVLRIGNNKEVQVSAFCLKRCQYLHEVELTTLNFMNVWKLSSSSHPGSEAPESNGLHRWWQDLCSVFA
TNDKLEVLMTNSVLGPPFLKALAAALRHPQCKLQKLLLRVNSTMLNQDLIGVLTGNQHLYLEIQHVE
VESKAVKLLCRVLRSPRCRLQCLRLEDCLATPRIWTDLGNNLQGNHGLKTLILRKNLENCGAYYLSVAQ
LERLSIENCNLTQLTCESLASCLRQSKMLTHLSLAENALKDEGAKHIWNALPHLRCPQLRVLVRKCDLTF
NCCQDMISALCKNKTLSLDLSFNLSKDDGVILLCEALKNPDCTLQILELENCLFTSICCCQAMASMLRKN
QHRLHLDLSKNAIGVYGILTLCFAFSSQKKREEVIFCIPAWTRITFSPTPHPPDFTGKSDCLSQINP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag: C-Myc/DDK
Predicted MW: 119.2 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_789781
Locus ID:	126205
UniProt ID:	Q86W28
RefSeq Size:	3934
Cytogenetics:	19q13.43
RefSeq ORF:	3144
Synonyms:	CLR19.2; NALP8; NOD16; PAN4
Summary:	This gene encodes a member of the nucleotide-binding oligomerization domain/ leucine rich repeat/ pyrin domain containing (NLRP) subfamily, which belongs to the Nod-like receptor family of proteins. NLRP genes play roles in the mammalian innate immune system through inflammasome formation and activation of caspases. In addition, NLRP genes have been found to function during mammalian reproduction. Consistent with a function during human preimplantation development, this gene is expressed at high levels in oocytes with decreased levels in embryos. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]

Product images:



Coomassie blue staining of purified NLRP8 protein (Cat# [TP318775]). The protein was produced from HEK293T cells transfected with NLRP8 cDNA clone (Cat# [RC218775]) using MegaTran 2.0 (Cat# [TT210002]).