

Product datasheet for TP318781L

PMP22 (NM_153322) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human peripheral myelin protein 22 (PMP22), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218781 representing NM_153322 Red =Cloning site Green =Tags(s)
	 MLLLLLSIIVLHVAVLVLLFVSTIVSQWIVGNHATDLWQNCSTSSSGNVVHHCFS SSPNEWLQSVQATMI LSIIFSILSLFFCQLFTLTGGGRFYITGIFQILAGLCVMSAAAITYVRHPEWHLNSDYSYGFAYILAW VAFPLALLSGVIYVILRKRE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_696997
Locus ID:	5376
UniProt ID:	Q01453 , Q6FH25



[View online »](#)

RefSeq Size: 1669

Cytogenetics: 17p12


RefSeq ORF: 480

Synonyms: CIDP; CMT1A; CMT1E; DSS; GAS-3; GAS3; HMSNIA; HNPP; Sp110

Summary: This gene encodes an integral membrane protein that is a major component of myelin in the peripheral nervous system. Studies suggest two alternately used promoters drive tissue-specific expression. Various mutations of this gene are causes of Charcot-Marie-Tooth disease Type IA, Dejerine-Sottas syndrome, and hereditary neuropathy with liability to pressure palsies. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified PMP22 protein (Cat# [TP318781]). The protein was produced from HEK293T cells transfected with PMP22 cDNA clone (Cat# [RC218781]) using MegaTran 2.0 (Cat# [TT210002]).