

Product datasheet for **TP318952L**

MAGEA4 (NM_001011548) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens melanoma antigen family A, 4 (MAGEA4), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218952 representing NM_001011548 Red=Cloning site Green=Tags(s)

MSSEQKSQHCKPEEGVEAQEEALGLVGAQAPTTEEQEA AVSSSSPLVPGTLEEVPAAESAGPPQSPQGAS
ALPTTISFTCWQPNEGSSSQEEEGPSTSPDAESLFREALSNKVDELAHFLLRKYRAKELVTKAEMLERV
IKNYKRCFPVIFGKASESLKMIFGIDVKEVDPASNTYTLVTCGLGSDGLLGNNQIFPKTGLLIIVLGTI
AMEGDSASEEEIWEELGVMGVYDGREHTVYGEPRKLLTQDWVQENYLEYRQVPGSNPARYEFLWGPRALA
ETSYVKVLEHVWRVNRVRIAYPSLREAALLEEEEGV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	34.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	ELISA capture for autoantibodies (PMID: 27323861) ELISA capture for autoantibodies (PMID: 27793776)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.




[View online »](#)

RefSeq:	NP_001011548
Locus ID:	4103
UniProt ID:	P43358 , A0A024RC12
RefSeq Size:	1724
Cytogenetics:	Xq28
RefSeq ORF:	951
Synonyms:	CT1.4; MAGE-41; MAGE-X2; MAGE4; MAGE4A; MAGE4B

Summary: This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. Several variants encoding the same protein have been found for this gene. [provided by RefSeq, Aug 2020]

Product images:



Coomassie blue staining of purified MAGEA4 protein (Cat# [TP318952]). The protein was produced from HEK293T cells transfected with MAGEA4 cDNA clone (Cat# [RC218952]) using MegaTran 2.0 (Cat# [TT210002]).