

Product datasheet for TP318977

UBE2D3 (NM_181893) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitin-conjugating enzyme E2D 3 (UBC4/5 homolog, yeast) (UBE2D3), transcript variant 9, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC218977 representing NM_181893 Red =Cloning site Green =Tags(s)
	<p>MLSNRKCLSKELSDLARDPPAQCSAGPVGDDMFHWQATIMGPNDSPLYQGGVFFLTIHFPTDYPFKPPKV A FTTRIYHPNINSNGSICLDILRSQWSPALTISKVLLSICLLCDPNPDDPLVPEIARIYKTDRDKYNRIS REWTQKYAM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	16.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_871622</u>
Locus ID:	7323




[View online »](#)

UniProt ID: [P61077](#)
RefSeq Size: 2006
Cytogenetics: 4q24
RefSeq ORF: 447
Synonyms: E2(17)KB3; UBC4/5; UBCH5C

Summary: The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme functions in the ubiquitination of the tumor-suppressor protein p53, which is induced by an E3 ubiquitin-protein ligase. [provided by RefSeq, Jan 2017]

Protein Pathways: Ubiquitin mediated proteolysis

Product images:



Coomassie blue staining of purified UBE2D3 protein (Cat# TP318977). The protein was produced from HEK293T cells transfected with UBE2D3 cDNA clone (Cat# [RC218977]) using MegaTran 2.0 (Cat# [TT210002]).