

Product datasheet for **TP319152M**

CAMKK2 (NM_153499) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calcium/calmodulin-dependent protein kinase kinase 2, beta (CAMKK2), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219152 representing NM_153499 Red =Cloning site Green =Tags(s)

MSSCVSSQPSSNRAAPQDELGGRGSSSESQKPCEALRGLSSLSIHLGMESFIVVTECEPGCAVDLGLAR
DRPLEADGQEVPLDSSGSQARPHLSGRKLSLQERSQGGLAAGGSLDMNGRCICPSLPYSPVSSPQSSPRL
PRRPTVESHHSITGMQDCVQLNQYTLKDEIGKGSYGVVKLAYNENDNTYYAMKVLSSKKLIRQAGFPRR
PPPRGTRPAPGGCIQPRGPIEQVYQEIAILKCLDHPNVVVKLVEVLDDPNEDHLYMVFELVNQGPVMEVPT
LKPLSEDQARFYFQDLIKGIEYLHYQKIIHRDIKPSNLLVGEDGHIKIADFGVSNEFKGSDALLSNTVGT
PAFMAPESLSETRKIFSGKALDVWAMGVTLYCFVFGQCPFMDERIMCLHSKIKSQALEFPDQPDIAEDLK
DLITRMLDKNPESRIVVPEIKLHPWVTRHGAEPLSEDENCTLVEVTEEEVENSVKHIPSLATVILVKTM
IRKRSFGNPFEGSRREERSLSAPGNLLTKQGSSENLQGTDPVVGGEEVLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	59.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_705719
Locus ID:	10645
UniProt ID:	Q96RR4 , A0A024RBP6
RefSeq Size:	5577
Cytogenetics:	12q24.31
RefSeq ORF:	1623
Synonyms:	CAMKK; CAMKKB
Summary:	The product of this gene belongs to the Serine/Threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. The major isoform of this gene plays a role in the calcium/calmodulin-dependent (CaM) kinase cascade by phosphorylating the downstream kinases CaMK1 and CaMK4. Protein products of this gene also phosphorylate AMP-activated protein kinase (AMPK). This gene has its strongest expression in the brain and influences signalling cascades involved with learning and memory, neuronal differentiation and migration, neurite outgrowth, and synapse formation. Alternative splicing results in multiple transcript variants encoding distinct isoforms. The identified isoforms differ in their ability to undergo autophosphorylation and to phosphorylate downstream kinases. [provided by RefSeq, Jul 2012]
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	Adipocytokine signaling pathway

Product images:



Coomassie blue staining of purified CAMKK2 protein (Cat# [TP319152]). The protein was produced from HEK293T cells transfected with CAMKK2 cDNA clone (Cat# [RC219152]) using MegaTran 2.0 (Cat# [TT210002]).