

Product datasheet for TP319337

CCDC46 (CEP112) (NM_001037325) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 46 (CCDC46), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219337 representing NM_001037325 Red =Cloning site Green =Tags(s)

MWASLSLDHPSAKENQALRLIEMREENGVPKTEQAGSLKPLRDTGKSNLKEKKANSKCLKQIEKEYTQKL
AKSSQIIAELQTTISSLKEENSQQQLAAERRLQDVRQKFEDEKKQLIRDNDQAIKVLQDELENRSNQVRC
AEKKLQHKELSESQEQITYIRQEYETKLKGLMPASLRQELEDTISSLSQVNFQKRASILQEELTTYQGR
R

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	24.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001032402
Locus ID:	201134




[View online »](#)

UniProt ID: [Q8N8E3](#)
RefSeq Size: 1258
Cytogenetics: 17q24.1
RefSeq ORF: 633
Synonyms: CCDC46; MACOCO; SPGF44

Summary: This gene encodes a coiled-coil domain containing protein that belongs to the cell division control protein 42 effector protein family. In neurons, it localizes to the cytoplasm of dendrites and is also enriched in the nucleus where it interacts with the RNA polymerase III transcriptional repressor Maf1 to regulate gamma-aminobutyric acid A receptor surface expression. In addition, the protein has been identified as a component of the human centrosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2014]

Product images:



Coomassie blue staining of purified CEP112 protein (Cat# TP319337). The protein was produced from HEK293T cells transfected with CEP112 cDNA clone (Cat# [RC219337]) using MegaTran 2.0 (Cat# [TT210002]).