

Product datasheet for **TP319357L**

AGXT2 (NM_031900) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human alanine-glyoxylate aminotransferase 2 (AGXT2), nuclear gene encoding mitochondrial protein, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219357 representing NM_031900 Red =Cloning site Green =Tags(s)

MTLIWRHLLRPLCLVTSAPRILEMHPFLSLGTSRTSVTKLSLHTKPRMPPCDFMPERYQSLGYNRVLEIH
KEHLSPVVTAYFQKPLLLHQGHMEWLFDAEGSRYLDFFSGIVTVSVGHCHPKVNAVAQKQLGRLWHTSTV
FFHPPMHEYAEKLAALLPEPLKVIFLVNSGSEANELAML MARAHSNNIDIISFRGAYHGCSPYTLGLTNV
GIYKMELPGGTGCQPTMCPDVFRGPWGGSHCRDSPVQTIRKCSAPDCCQAKDQYIEQFKDTLSTVAKS
IAGFFAEPIQGVNGVWQYPKGFLEAFELVRARGGVCIADDEVQTFGRLGSHFWGFQTHDVLDPDIVTMAK
GIGNGFPMMAAVITTP EIAKSLAKCLQHFNTFGGNP MACAIGSAVLEVIKEENLQENSQEVGTYMLLKFAK
LRDEFEIVGDVRGKGLMIGIEMVQDKISCRPLPREEVNIHEDCKHMGLLVGRGSIFSQTFRIAPSMCIT
KPEVDFAVEVFRSALTQHMERRAK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	52.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_114106
Locus ID:	64902
UniProt ID:	Q9BYV1
RefSeq Size:	2165
Cytogenetics:	5p13.2
RefSeq ORF:	1542
Synonyms:	AGT2; BAIBA; DAIBAT
Summary:	The protein encoded by this gene is a class III pyridoxal-phosphate-dependent mitochondrial aminotransferase. It catalyzes the conversion of glyoxylate to glycine using L-alanine as the amino donor. It is an important regulator of methylarginines and is involved in the control of blood pressure in kidney. Polymorphisms in this gene affect methylarginine and beta-aminoisobutyrate metabolism, and are associated with carotid atherosclerosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism

Product images:



Coomassie blue staining of purified AGXT2 protein (Cat# [TP319357]). The protein was produced from HEK293T cells transfected with AGXT2 cDNA clone (Cat# [RC219357]) using MegaTran 2.0 (Cat# [TT210002]).