

Product datasheet for TP319368L

NMDAR1 (GRIN1) (NM_000832) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutamate receptor, ionotropic, N-methyl D-aspartate 1 (GRIN1), transcript variant NR1-1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219368 representing NM_000832 Red=Cloning site Green=Tags(s)

MSTMRLTLALLFSCSVARAACDPKIVNIGAVLSTRKHEQMFREAVNQANKRHGSWKIQLNATSVTHKPN
AIQMALSVCEDLISSQVYAILVSHPTPNDFHTPTVSYTAGFYRIPVLGLTTRMSIYSDKSIHLSFLRT
VPPYSHQSSVWFEMMRVYSWNHIILLVSDDHEGRAAQKRLETLEERESKAEKVLQFDPGKNTVALLME
AKELEARVIILSASEDDAATVYRAAAMLNMTGSGYVWLVEREISGNALRYAPDGILGLQLINGKNESAH
ISDAVGVAQAVHELLEKENITDPPRGCVGNTNIWKTGPLFKRVLMSKYADGVTGRVEFNEDGDRKFAN
YSIMNLQNRKLVQVGIYNGTHVIPNDRKIIWPGGETEKPRGYQMSTRLKIVTIHQEPFVYVKPTLSDGTC
KEEFTVNGDPVKKVICTGPNDTSPGSPRHTVPQCCYGFICIDLLIKLARTMNFTEVHLVADGKFGTQERV
NNSNKKEWNGMMGELLSGQADMIVAPLTINNERAQYIEFSKPFKYQGLTILVKEIPRSTLDSFMQPFQS
TLWLLVGLSVHVAVMLYLLDRFSPFGRFKVNSEEEEEEDALTSSAMWFSWGVLNLSGIGEGAPRSFSAR
ILGMVWAGFAMIIVASYTANLAAFLVDRPEERITGINDPRLRNPSPDKFIYATVKQSSVDIYFRRQVELS
TMYRHMEKHNYESAAEAIQAVRDNKLHAFIWDASAVLEFEASQKCDLVTTGELFFRSFGFYGMRKDSWPWKQ
NVLSILKSHENGFMEDLDKTWVRYQECDSRSNAPATLTFENMAGVFMVAGGIVAGIFLIFIEIAYKRH
KDARRKQMLAFAAVNVWRKNLQQYHPTDITGPLNLSDPVSTVV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	97.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000823
Locus ID:	2902
UniProt ID:	Q05586
RefSeq Size:	3902
Cytogenetics:	9q34.3
RefSeq ORF:	2655
Synonyms:	GluN1; MRD8; NDHMSD; NDHMSR; NMD-R1; NMDA1; NMDAR1; NR1
Summary:	The protein encoded by this gene is a critical subunit of N-methyl-D-aspartate receptors, members of the glutamate receptor channel superfamily which are heteromeric protein complexes with multiple subunits arranged to form a ligand-gated ion channel. These subunits play a key role in the plasticity of synapses, which is believed to underlie memory and learning. Cell-specific factors are thought to control expression of different isoforms, possibly contributing to the functional diversity of the subunits. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane
Protein Pathways:	Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Calcium signaling pathway, Huntington's disease, Long-term potentiation, Neuroactive ligand-receptor interaction

Product images:



Coomassie blue staining of purified GRIN1 protein (Cat# [TP319368]). The protein was produced from HEK293T cells transfected with GRIN1 cDNA clone (Cat# [RC219368]) using MegaTran 2.0 (Cat# [TT210002]).