

## Product datasheet for **TP319375L**

### LCK (NM\_005356) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human lymphocyte-specific protein tyrosine kinase (LCK), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	Recombinant protein was produced with TrueORF clone, RC219375.
Tag:	C-Myc/DDK
Predicted MW:	57.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_005347</a>
Locus ID:	3932
UniProt ID:	<a href="#">P06239</a> , <a href="#">A0A0S2Z3Y8</a>
RefSeq Size:	2032
Cytogenetics:	1p35.2
RefSeq ORF:	1527
Synonyms:	IMD22; LSK; p56lck; pp58lck; YT16



[View online »](#)

**Summary:**


This gene is a member of the Src family of protein tyrosine kinases (PTKs). The encoded protein is a key signaling molecule in the selection and maturation of developing T-cells. It contains N-terminal sites for myristylation and palmylation, a PTK domain, and SH2 and SH3 domains which are involved in mediating protein-protein interactions with phosphotyrosine-containing and proline-rich motifs, respectively. The protein localizes to the plasma membrane and pericentrosomal vesicles, and binds to cell surface receptors, including CD4 and CD8, and other signaling molecules. Multiple alternatively spliced variants encoding different isoforms have been described. [provided by RefSeq, Aug 2016]

**Protein Families:**

Druggable Genome, Protein Kinase, Stem cell - Pluripotency

**Protein Pathways:**

Natural killer cell mediated cytotoxicity, Primary immunodeficiency, T cell receptor signaling pathway

**Product images:**

Coomassie blue staining of purified LCK protein (Cat# [TP319375]). The protein was produced from HEK293T cells transfected with LCK cDNA clone (Cat# [RC219375]) using MegaTran 2.0 (Cat# [TT210002]).