

Product datasheet for **TP319393**

Ninjurin2 (NINJ2) (NM_016533) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ninjurin 2 (NINJ2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219393 representing NM_016533 Red =Cloning site Green =Tags(s)
	 MAGLSRQLCALSHPKKAAETQTAEPGGAHAVCSRHPVRVKGLEGSEMESARENIDLQPGSSDPRSQPINL NHYATKKSVAESMLDVALFMSNAMRLKAVLEQGPSSHYTTLVTLISLSSLLQWIGVLLVVIARLNLNE VEKQWRLNQLNNAATILVFFTWINVFITAFGAHKTGFLAARASNPL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057617
Locus ID:	4815
UniProt ID:	Q9NZG7 , A0A590UJR9 , B4DJC1
RefSeq Size:	1073



[View online »](#)


Cytogenetics: 12p13.33

RefSeq ORF: 564

Summary: The protein encoded by this gene belongs to the ninjurin (for nerve injury induced) family. It is a cell surface adhesion protein that is upregulated in Schwann cells surrounding the distal segment of injured nerve, and promotes neurite outgrowth, thus may have a role in nerve regeneration after nerve injury. [provided by RefSeq, Oct 2011]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified NINJ2 protein (Cat# TP319393). The protein was produced from HEK293T cells transfected with NINJ2 cDNA clone (Cat# [RC219393]) using MegaTran 2.0 (Cat# [TT210002]).