

Product datasheet for TP319422L

HVCN1 (NM_032369) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydrogen voltage-gated channel 1 (HVCN1), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219422 protein sequence Red=Cloning site Green=Tags(s)

MATWDEKAVTRRAKVAPAERMSKFLRHFTVGGDDYHAWNINYYKKWENEEEEEEEEQPPPTPVS GEEGRAA
APDVAPAPGPAPRAPLDFRGMRLKLFSSHRFQVIIIICLVVLDALLVLAELIDLKIIQPKNNYAAMVFH
YMSITILVFFMMEIIFKLFVFRLEFFHHKFEILDVVVVVSFILDIVLLFQEHQFEALGLLILLRLWRVA
RIINGIIISVKTRSERQLLRKQMNVQLAAKIQHLEFSCSEKEQEIERLNKLLRQHGLLGEVN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	31.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115745
Locus ID:	84329



[View online »](#)

UniProt ID: [Q96D96](#)

RefSeq Size: 1727

Cytogenetics: 12q24.11


RefSeq ORF: 819

Synonyms: HV1; VSOP

Summary: This gene encodes a voltage-gated protein channel protein expressed more highly in certain cells of the immune system. Phagocytic cells produce superoxide anions which require this channel protein, and in B cells this same process facilitates antibody production. This same channel protein, however, can also regulate functions in other cells including spermatozoa. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]

Protein Families: Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified HVCN1 protein (Cat# [TP319422]). The protein was produced from HEK293T cells transfected with HVCN1 cDNA clone (Cat# [RC219422]) using MegaTran 2.0 (Cat# [TT210002]).