

Product datasheet for **TP319432**

BAGE2 (NM_182482) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human B melanoma antigen family, member 2 (BAGE2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219432 representing NM_182482 Red =Cloning site Green =Tags(s)
	MAAGVVFLALSAQLLQARLMKEESPVSWRLEPEDGTALDVHVFVSTLEPLSNVAVKRNVPKCIILVLQEP TPFRISVTSSCFVQNTLTLLKDRRKMQTVQCATARETS
	TR TRPLE QKLISEEDLA NDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	12 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_872288
Locus ID:	85319
RefSeq Size:	1891
Cytogenetics:	21p11.2
RefSeq ORF:	327




[View online »](#)

Synonyms: CT2.2

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified BAGE2 protein (Cat# TP319432). The protein was produced from HEK293T cells transfected with BAGE2 cDNA clone (Cat# [RC219432]) using MegaTran 2.0 (Cat# [TT210002]).