

Product datasheet for **TP319442L**

FGD4 (NM_139241) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human FYVE, RhoGEF and PH domain containing 4 (FGD4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219442 representing NM_139241 Red =Cloning site Green =Tags(s)

MEEIKPASASCVSKEKPSKVSDLISRFEGGSSLSNYSDLKKESAVNLNAPRTPGRHGLTTTPQQKLLSQH
LPQRQGNDDTKTQGAQTCVANGVMAAQNQMECEEEKAATLSSDTSIQASEPLLDMHIVNGERDETATAPA
SPTTDSCDGNASDSSYRTPGIGPVLPLEERGAETETKVQERENGESPLELEQLDQHHMKETNEQKLHKI
ANELLTERAYVNRLLDLDQVFYCKLLEEANRGSFPAEMVNKIFSNISSINAFHSKFLLEPELEKRMQWE
TTPRIGDILQKLAPFLKMYGEYVKGFDNAMELVKNMTERIPQFKSVVEIQKQKICGSLTLQHMHLEPVQ
RIPRYEMLLKDYLRLKPPDSLWDNDAKKSLEIISTAASHSNSAIRKMNKLLLEIYEMLGEEEDIVNPS
NELIKEGQILKLAARNTSAQERYLFLFNNMLLYCVPKFSLVGSKFTVRTRVIGIDGMKIVETQNEEYPHTF
QVSGKERTLELQASSAQDKEEWIKALQETIDAFHQRHETFRNAIAKDNDIHSEVSTAELGKRAPRWIRDN
EVTMCMKCKEPFNALTRRRHHCACGYVVCWKCSYKAQLEYDGGKLSKVCKDCYQIISGFTDSEEKRRK
GILEIESAEVSGNSVCSFLQYMEKSKPWQKAWCVIPKQDPLVLYMYGAPQDVRAQATIPLLGYVDEMP
RSADLPHSFKLTQSKSVHSFAADSEELKQKWLKVILLAVTGETPGGPNEHPATLDDHPEPKKKSEC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	86.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_640334
Locus ID:	121512
UniProt ID:	Q96M96
RefSeq Size:	2931
Cytogenetics:	12p11.21
RefSeq ORF:	2298
Synonyms:	CMT4H; FRABP; ZFYVE6
Summary:	This gene encodes a protein that is involved in the regulation of the actin cytoskeleton and cell shape. This protein contains an actin filament-binding domain, which together with its Dbl homology domain and one of its pleckstrin homology domains, can form microspikes. This protein can activate MAPK8 independently of the actin filament-binding domain, and it is also involved in the activation of CDC42 via the exchange of bound GDP for free GTP. The activation of CDC42 also enables this protein to play a role in mediating the cellular invasion of <i>Cryptosporidium parvum</i> , an intracellular parasite that infects the gastrointestinal tract. Mutations in this gene can cause Charcot-Marie-Tooth disease type 4H (CMT4H), a disorder of the peripheral nervous system. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2015]

Product images:



Coomassie blue staining of purified FGD4 protein (Cat# [TP319442]). The protein was produced from HEK293T cells transfected with FGD4 cDNA clone (Cat# [RC219442]) using MegaTran 2.0 (Cat# [TT210002]).