

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_060628](#)

Locus ID: 22950

UniProt ID: [Q9BWU0](#)


RefSeq Size: 2982

Cytogenetics: 2p23.3

RefSeq ORF: 2388

Synonyms: HLC3

Product images:



Coomassie blue staining of purified SLC4A1AP protein (Cat# [TP319504]). The protein was produced from HEK293T cells transfected with SLC4A1AP cDNA clone (Cat# [RC219504]) using MegaTran 2.0 (Cat# [TT210002]).