

Product datasheet for TP319561

H2BC17 (NM_003527) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human histone cluster 1, H2bo (HIST1H2BO), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219561 protein sequence Red =Cloning site Green =Tags(s)
	MPDPAKSAPAPKKGSKKAVTKAQKKDGGKRRKRSRKESYSIYVYKVLKQVHPDTGISSKAMGIMNSFVNDI FERIAGEASRLAHYNKRSTITSREIQTAVRLLLPGELAKHAVSEGTKAVTKYTSSK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003518
Locus ID:	8348
UniProt ID:	P23527
RefSeq Size:	467
Cytogenetics:	6p22.1



[View online »](#)

RefSeq ORF: 378


Synonyms: dj193B12.2; H2B.2; H2B/n; H2BFN; HIST1H2BO

Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2B family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the small histone gene cluster on chromosome 6p22-p21.3. [provided by RefSeq, Aug 2015]

Protein Families: Druggable Genome

Protein Pathways: Systemic lupus erythematosus

Product images:



Coomassie blue staining of purified HIST1H2BO protein (Cat# TP319561). The protein was produced from HEK293T cells transfected with HIST1H2BO cDNA clone (Cat# [RC219561]) using MegaTran 2.0 (Cat# [TT210002]).