

Product datasheet for **TP319593L**

DGKH (NM_178009) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human diacylglycerol kinase, eta (DGKH), transcript variant 2, 1 mg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC219593 representing NM_178009 Red =Cloning site Green =Tags(s) |

MAGAGGQHPPGAAGGAAAGAGAAVTSAAASAGPGEDSSDSEAEQEGPQKLIRKVSTSGQIRTKTSIKEG
QLLKQTSSFQRWKKRYFKLRGRTLYYAKDSKSLIFDEVDLSDASVAEASTKNANNSTIITPFRRMLLCA
ENRKEMEDWISSLKSVQTPREPYEVAQFNVEHFGSMHNWYACSHARPTFCNVCRESLSGVTSHGLSCEVCK
FKAHKRAVRATNNCKWTTLASIGKDIEDEDGVAMPHQWLEGNLPVSAKCAVCDKTCGVLRLQDWKCL
WCKTMVHTACKDLYHPICPLGQCKVSIIPPIALNSTSDGFCRATFSFCVSPLLVFNKSGDNQGVKFL
RRFKQLLNPAQVFDLMNGGPHLGLRFLQKFDNFRILVCGGDGSGVWVWVSEIDKLNLNKQCQLGVPLGTG
NDLARVLGWGGSYDDDTQLPQILEKLERASTKMLDRWSIMTYELKLPKASLLPGPPEASEEFYMTIYED
SVATHLTKILNSDEHAVVISSAKTLCETVKDFVAKVEKTYDKTLENVAVADAVASKCSVLNEKLEQLLQA
LHTDSQAAPVLPGLSPLIVEEDAVESSSEESLGESEKQLGDDVTKPSSQKAVKPREIMLRANSLKKAQRQ
VIEEAGKVMDDPTVHPCEPANQSSDYDSTETDESKEEAKDDGAKESITVKTAPRSPDARASYGHSQTDVSV
PGPAVAASKENLPVLNTRICPGLRAGLAASIAGSSIINKMLLANIDPFGATPFIDPDLDSDVDGYSEKCV
MNNYFGIGLDAKISLEFNKREEHPEKCRSRTKNLMWYGVLTRELLQRSYKNLEQRVQLECDGQYIPLP
SLQGIAVLNIPSYAGGTNFWGGTKEDDIFAAPSFDDKILEVVAIFDSMQMAVSRVIKLVHRIAQCRTVK
ITIFGDEGVPVQVDGEAWVQPPGIIKIVHKNRAQMLTRDRAFESTLKSWEKQKCDGSKPVLRTHLYIHH
AIDLATEEVSQMQLCSQAAEELITRICDAATIHCLLEQELAHAVNACSHALNKANPRCPESLTRDTATEI
AINVKALYNETESLLVGRVPLQLESPHEERVSNALHSVEVELQKLTEIPWLYIILHPNEDEEPPMDCTKR
NNRSTVFRIVPKFKKEKVQKQKTSSQPVQKWGTEEVAAWLDLLNLGEYKDFIRHDIRGAELLHLERRDL
KDLGIPKVGHVKRILQGIKELGRSTPQSEV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


| | |
|----------------|---|
| Tag: | C-Myc/DDK |
| Predicted MW: | 134.7 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |



[View online »](#)

| | |
|--------------------------|---|
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_821077 |
| Locus ID: | 160851 |
| UniProt ID: | Q86XP1 |
| RefSeq Size: | 4248 |
| Cytogenetics: | 13q14.11 |
| RefSeq ORF: | 3660 |
| Synonyms: | DGKeta |
| Summary: | This gene encodes a member of the diacylglycerol kinase (DGK) enzyme family. Members of this family are involved in regulating intracellular concentrations of diacylglycerol and phosphatidic acid. Variation in this gene has been associated with bipolar disorder. Alternatively spliced transcript variants have been identified. [provided by RefSeq, Jul 2014] |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Phosphatidylinositol signaling system |

Product images:



Coomassie blue staining of purified DGKH protein (Cat# [TP319593]). The protein was produced from HEK293T cells transfected with DGKH cDNA clone (Cat# [RC219593]) using MegaTran 2.0 (Cat# [TT210002]).