

Product datasheet for **TP319625L**

COX8C (NM_182971) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cytochrome c oxidase subunit 8C (COX8C), nuclear gene encoding mitochondrial protein, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219625 protein sequence Red =Cloning site Green =Tags(s)
	MPLLRGRCPARRHYRRLALLGLQPAPRFAHSGPPRQRPLSAAEMAVGLVVFVFTTFLTPAAYVLGNLKQFR RN
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	7.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_892016
Locus ID:	341947
UniProt ID:	Q7Z4L0
RefSeq Size:	531



[View online »](#)

Cytogenetics:	14q32.12
RefSeq ORF:	216
Synonyms:	COX8-3
Summary:	This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transmembrane
Protein Pathways:	Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images:



Coomassie blue staining of purified COX8C protein (Cat# [TP319625]). The protein was produced from HEK293T cells transfected with COX8C cDNA clone (Cat# [RC219625]) using MegaTran 2.0 (Cat# [TT210002]).