

Product datasheet for **TP319637L**

CPNE8 (NM_153634) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human copine VIII (CPNE8), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219637 protein sequence Red =Cloning site Green =Tags(s)

MDSRYNSTAGIGDLNQLSAAIPATRVEVSVSCRNLLDRDTFSKSDPICVLYVQGVGNKEWREFGRTEVID
NTLNPDFVRKFILDYFFEERENLRFDLVDVDSKSPNLSKHDFLGQVFCTLGEIVGSQGSRLKPIVGIPG
KKCGTIILTAELNCCRDVLMQFCANKLDDKDFGKSDPFLVFYRSNEDGSFTICHKTEVVKNTLNPVW
QAFKISVRALCNGDYDRTIKVEYDWD RDGSHDFIGFTTSYRELSRGQSQFNVEVWNP KKKGKKKKYT
NSGTVTLLSFLVETEVSFLDYIKGGTQINFTVAIDFTASNGNPAQPTSLHYMNPYQLNAYGMALKAVGEI
VQDYDSDKMFPALGFGAKLPDGRISHEFALNGNPQNPYCDGIEGVMEAYYRSLKSVQLYGPTNFAPVIN
HVARYASSVKDGSQYFVLLIVTDGVISDMAQTKEIVNASKLPMSIIVGVGPAEFDAMVELDGD D V R V S
SRGKYAERDQVFPFRDYIDRSNGHILSMARLAKDVLAIEPQFLSYMRARGIKPSPAPPPYTPPTHVL
QTQI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	62.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_705898</u>
Locus ID:	144402
UniProt ID:	<u>Q86YQ8</u>
RefSeq Size:	3476
Cytogenetics:	12q12
RefSeq ORF:	1692
Summary:	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene is one of several genes that encode a calcium-dependent protein containing two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified CPNE8 protein (Cat# [TP319637]). The protein was produced from HEK293T cells transfected with CPNE8 cDNA clone (Cat# [RC219637]) using MegaTran 2.0 (Cat# [TT210002]).