

Product datasheet for TP319680

CCDC159 (NM_001080503) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain-containing-like (LOC126075), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219680 representing NM_001080503 Red =Cloning site Green =Tags(s)
	 MGEHEQVKPLETSSSKVKAKTIVMIPDSQKLLRCELESLKSQQLQAQTKAFEFLNHSVTMLEKESCLQQIK IQQLEEVLSPTGRQGEKEEHKWGMEQGRQELYGALTQGLGLEKTLRDSEEMQRARTTRCLQLLAQEIRD SKKFLWEELELVREEVTFIYQKLQAQEQEISENLVNIQKMQKTQVKCRKILTKMKQQGHETAACPETEEI PQGASGCWKDDLQKELSDIWSAVHVLQNSIDSLTLCGACPKASSLRGHKGGHQLSPPLPSWDSDDSDCDQ DLSQPPFSKSGRSFPPA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	33.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001073972
Locus ID:	126075



[View online »](#)

UniProt ID: [POC716](#)
RefSeq Size: 2238
Cytogenetics: 19p13.2
RefSeq ORF: 891

Product images:



Coomassie blue staining of purified CCDC159 protein (Cat# TP319680). The protein was produced from HEK293T cells transfected with CCDC159 cDNA clone (Cat# [RC219680]) using MegaTran 2.0 (Cat# [TT210002]).