

Product datasheet for **TP319692**

TPD52L2 (NM_199359) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human tumor protein D52-like 2 (TPD52L2), transcript variant 6, 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC219692 representing NM_199359
Red=Cloning site **Green**=Tags(s)

MDSAGQDINLNNSPNKGLLSDSMTDVPVDTGVAARTPAVEGLTEAAAAEELRAELTKVEEIVTLRQVLA
AKERHCGELKRRRLGLSTLGELKQNLRSRSHDVQVSSAYKKTQETLSQAGQKTSALSTVGSASIRKLGDMRN
SATFKSFEDRVGTIKSKVVGDRENGSDNLPSSAGSGDKPLSDPAPF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 19.7 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_955391](#)
Locus ID: 7165
UniProt ID: [Q68E05](#)
RefSeq Size: 2297



[View online »](#)


Cytogenetics: 20q13.33

RefSeq ORF: 559

Synonyms: D54; TPD54

Summary: This gene encodes a member of the tumor protein D52-like family. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. Expression of this gene may be a marker for breast cancer and acute lymphoblastic leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 12. [provided by RefSeq, Aug 2011]

Product images:



Coomassie blue staining of purified TPD52L2 protein (Cat# TP319692). The protein was produced from HEK293T cells transfected with TPD52L2 cDNA clone (Cat# [RC219692]) using MegaTran 2.0 (Cat# [TT210002]).