

Product datasheet for **TP319746L**

SFRS12IP1 (SREK1IP1) (NM_173829) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human SFRS12-interacting protein 1 (SFRS12IP1), 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or >RC219746 protein sequence

AA Sequence: **Red**=Cloning site **Green**=Tags(s)

MAVPGCNKDSVRAGCKKCGYPGHILTFECRNFLRVDPKRDIVLDVSSTSSSEDSDEENEELNKLQALQEKRI
NEEEKKKEKSKEKIKLKKKRKRSYSSSSTEEDTSKQKKQKYQKKEKKKEKSSKSKGKHHKKEKKRKK
EKHSSTPNSSEFSRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 18 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_776190](#)



Locus ID:	285672
UniProt ID:	<u>Q8N9Q2</u>
RefSeq Size:	6856
Cytogenetics:	5q12.3
RefSeq ORF:	465
Synonyms:	PI8SRP; SFRS12IP1
Summary:	Possible splicing regulator involved in the control of cellular survival.[UniProtKB/Swiss-Prot Function]

Product images:


Coomassie blue staining of purified SREK1IP1 protein (Cat# [TP319746]). The protein was produced from HEK293T cells transfected with SREK1IP1 cDNA clone (Cat# [RC219746]) using MegaTran 2.0 (Cat# [TT210002]).