

## Product datasheet for TP319762

### FHL2 (NM\_001039492) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens four and a half LIM domains 2 (FHL2), transcript variant 5, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219762 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MTERFDCHHCNESLFGKKYILREESPYCVCFETLFANTCEECKPIGCDCKDLSYKDRHWHEACFHCSQ CRNSLVDKPFPAKEDQLLCTDCYSNEYSSKCQECKKTIMPGTRKMEYKGSSWHETCFICHRCQQPIGTKS FIPKDNQNFVCPCYEKQHAMQCVQCKPITGGVTYREQPWHKECFVCTACRKQLSGQRFTARDDFAYC L NCFCDLYAKKCAGCTNPISGLGGTKYISFEERQWHNDCFNCKKCSLSLVGRGFLTERDDILCPDCGKDI  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	32 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001034581</a></u>
Locus ID:	2274


[View online »](#)

UniProt ID:	<u>Q14192</u>
RefSeq Size:	1552
Cytogenetics:	2q12.2
RefSeq ORF:	837
Synonyms:	AAG11; DRAL; FHL-2; SLIM-3; SLIM3
Summary:	This gene encodes a member of the four-and-a-half-LIM-only protein family. Family members contain two highly conserved, tandemly arranged, zinc finger domains with four highly conserved cysteines binding a zinc atom in each zinc finger. This protein is thought to have a role in the assembly of extracellular membranes. Also, this gene is down-regulated during transformation of normal myoblasts to rhabdomyosarcoma cells and the encoded protein may function as a link between presenilin-2 and an intracellular signaling pathway. Multiple alternatively spliced variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2016]
Protein Families:	Druggable Genome

### Product images:



Coomassie blue staining of purified FHL2 protein (Cat# TP319762). The protein was produced from HEK293T cells transfected with FHL2 cDNA clone (Cat# [RC219762]) using MegaTran 2.0 (Cat# [TT210002]).