

Product datasheet for TP319857

PCDHGA12 (NM_003735) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protocadherin gamma subfamily A, 12 (PCDHGA12), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219857 representing NM_003735 Red=Cloning site Green=Tags(s)

MIPARLHRDYKGLVLLGILLGTLWETGCTQIRYSVPEELEKGSRVGDISRDLGLEPRELAERGVRIIPRG
RTQLFALNPRSGSLVTAGRIDREELCMGAIKCQLNLDILMEDKVKIYGVEVEVRDINDNAPYFRESELEI
KISENAATEMRFLPHAWDPDIGKNSLQSYELSPNTHFSLIVQNGADGSKYPELVLKRALDREEKAAHHL
VLTASDGGDPVRTGTARIRVMVLDANDNAPAFAPQPEYRASVPENLALGTQLLVNATDPDEGVNAEVRYS
FRYVDDKAAQVFKLDCNSGTISTIGELDHEESGFYQMEVQAMDNAGYSARAKVLITVLDVNDNAPEVLT
SLASSVPENSPRGTLIALLNVDQDSEENGQVICFIQGNLFPKLEKSYGNYSVLTDIVLDREQVPSYNI
TVTATDRGTPPLSTETHISLNVADTNDNPPVFPQASYSAYIPENNPRGVSLSVTAHDPDCEENAQITYS
LAENTIQQASLSSYVSINSDTGVLYALSSFDYEQFRDLQVKVMARDNGHPPLSSNVLSLFLVDQNDNAP
EILYPALPTDGGSTGVELAPRSAEPGYLVTKVAVDRDSGQNAWLSYRLLKASEPGLFSVGLHTGEVRTAR
ALLDRDALKQSLVAVQDHGQPPLSATVTLTVAVADSIPQVLADLGSLESPANSETSDLTYLVVAVAAV
SCVFLAFVILLALRLRRWHKSRLQASGGGLTGAPASHFVGDGVQAFQTYSHEVSLTTDSRKSHLIF
PQPNYADMLVSQESFEKSEPLLLSGDSVFSKDSHGLIEQAPPNTDWRFSQAQRPGTSGSQNGDDTGTWPN
NQFDTEMLQAMILASASEAADGSSTLGGGAGTMGLSARYGPQFTLQHVPDYRQNVYIPGSNATLTNAAGK
RDGKAPAGGNGNKKKSGKKEKK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	97.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003726
Locus ID:	26025
UniProt ID:	O60330
RefSeq Size:	4774
Cytogenetics:	5q31.3
RefSeq ORF:	2796
Synonyms:	CDH21; FIB3; PCDH-GAMMA-A12
Summary:	<p>This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]</p>
Protein Families:	Transmembrane

Product images:

Coomassie blue staining of purified PCDHGA12 protein (Cat# TP319857). The protein was produced from HEK293T cells transfected with PCDHGA12 cDNA clone (Cat# [RC219857]) using MegaTran 2.0 (Cat# [TT210002]).