

## Product datasheet for TP319929

### NVL (NM\_206840) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nuclear VCP-like (NVL), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219929 representing NM_206840 Red=Cloning site Green=Tags(s)

MEDYDPDQSANHMNSSLLSLYRKGPNPDSVSNTPEMEQRETTSSSTPRISSKTGSIPLKTPAKDSEGGWFID  
KTPSVKKDSFFLDLSCEKSNPKKPIEQDSKDSSLLESDMKRKGKLNKNGSKRKKEDLQEVGDGEIEAVL  
QKKAKARGLEFQISNVKFEDVGGNDMTLKEVCKMLIHMHRHPEVYHHLGVPPRGVLLHGPPGCGKTLAH  
AIAGELDLPILKVAAPEIVSGVSGESEQKLRELFEQAVSNAPCIIFIDEIDAITPKREVASKDMERRIVA  
QLLTCMDDLNNVAATARVLVIGATNRPDSDLPALRRAGRFDRICLGIPDEASRERILQTLCRKLRPQA  
FDFCHLAHLTPGFVGADLMALCREAAMCAVNRVLMKLQEQKKNPEMEDLPSKGVQEERLGTEPTSETQD  
ELQRLGLLRDQDPLSEEQMQGLCIELNDFIVALSSVQPSAKREGFVTPNVTWADIGALEDIREELTMA  
ILAPVRNPQFKALGLVTPAGVLLAGPPGCGKTLAKAVANESGLNFISVKPELLNMYVGESERAVRQV  
FQRAKNSAPCVIFFDEVDALCPRRSRETGASVRVQNQLLTEMDGLEARQQVFIMAATNRPDIIDPAILR  
PGRLDKTLFVGLPPPADRLAILKTITKNGTKPPLDADVNLEAIAGDLRCDCYTGADLSALVREASICALR  
QEMARQKSGNEKGELKVSHKHFEFAFKVRSSISKDKQIMYERLQESLSR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	82.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_996671</a>
<b>Locus ID:</b>	4931
<b>UniProt ID:</b>	<a href="#">O15381</a>
<b>RefSeq Size:</b>	2847
<b>Cytogenetics:</b>	1q42.11
<b>RefSeq ORF:</b>	2250
<b>Synonyms:</b>	NVL2
<b>Summary:</b>	This gene encodes a member of the AAA (ATPases associated with diverse cellular activities) superfamily. Multiple transcript variants encoding different isoforms have been found for this gene. Two encoded proteins, described as major and minor isoforms, have been localized to distinct regions of the nucleus. The largest encoded protein (major isoform) has been localized to the nucleolus and shown to participate in ribosome biosynthesis (PMID: 15469983, 16782053), while the minor isoform has been localized to the nucleoplasm. [provided by RefSeq, Aug 2011]

### Product images:



Coomassie blue staining of purified NVL protein (Cat# TP319929). The protein was produced from HEK293T cells transfected with NVL cDNA clone (Cat# [RC219929]) using MegaTran 2.0 (Cat# [TT210002]).