

Product datasheet for TP319963L

MZT2A (NM_001085365) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 128, member A (FAM128A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC219963 representing NM_001085365 Red=Cloning site Green=Tags(s)
	MELYELAQAAGGGIDPDVFKILVDLLKLNVAFLAVFQMLKSMCAGQRLASEPQDPAAVSLPTSSVPETRG RDKGSAAALGGVLALAERSNHEGSSQRMPRQPSATRLPKGGGPGKSPTQGST
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	16 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001078834
Locus ID:	653784
UniProt ID:	Q6P582
RefSeq Size:	708



[View online »](#)

Cytogenetics: 2q21.1
RefSeq ORF: 363
Synonyms: FAM128A; MOZART2A

Product images:

Coomassie blue staining of purified MZT2A protein (Cat# [TP319963]). The protein was produced from HEK293T cells transfected with MZT2A cDNA clone (Cat# [RC219963]) using MegaTran 2.0 (Cat# [TT210002]).