

Product datasheet for TP320009

FGFR1 (NM_023106) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fibroblast growth factor receptor 1 (FGFR1), transcript variant 4, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220009 protein sequence Red=Cloning site Green=Tags(s)

MWSWKCLLFWAVLVLTATLCTARPSPTLPEQDALPSEDDDDDDSSSEEKETDNTKPNPVAPYWTSPEKM
EKKLHAVPAAKTVKFKCPSSGTPNPTLRWLKNGKEFKPDHRIGGYKVRYATWSIIMDSVWPSDKGNYTCI
VENEYGSINHTYQLDVVERSPHRPILQAGLPANKTVALGSNVEFMCKVYSDPQPHIQWLKHIEVNGSKIG
PDNLPLYVQILKTAGVNTTDKEMEVLHLRNVSFEDAGEYTCLAGNSIGLSHSAWLTVLEALEERPVMVTS
PLYLEIIYCTGAFILISCMVGSVIVYKMKSGTKKSDFHQSMAVHKLAKSIPLRRQVTVSADSSASMNSGV
LLVRPSRLSSSGTPMLAGVSEYELPEDPRWELPRDRLVLGKPLGEGCFGQVLAEAIGLDKDKPNRVTKV
AVKMLKSDATEKDLSDLISEMEMMKMIGKHKNINLLGACTQDGPLYVIVEYASKGNLREYLQARRPPGL
EYCYNPSHNPEEQSSKDLVSCAYQVARGMEYLASKKCIHRDLAARNVLVTEDNVMKIADFLGARDIHHI
DYYKKTNGRLPVKWMPEALFDRIYTHQSDVWSFGVLLWEIFTLGGSPYGPVVEELFKLLKEGHRMDK
PSNCTNELYMMMRDCWHAVPSQRPTFKQLVEDLDRIVALTSNQEYLDLSMPLDQYSPSPDTRSSSTCSSG
EDSVFSHEPLPEEPCLPRHPAQLANGGLKRR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	79.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_075594
Locus ID:	2260
UniProt ID:	P11362
RefSeq Size:	5644
Cytogenetics:	8p11.23
RefSeq ORF:	2193
Synonyms:	bFGF-R-1; BFGFR; CD331; CEK; ECCL; FGFBR; FGFR-1; FLG; FLT-2; FLT2; HBGFR; HH2; HRTFDS; KAL2; N-SAM; OGD
Summary:	<p>The protein encoded by this gene is a member of the fibroblast growth factor receptor (FGFR) family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds both acidic and basic fibroblast growth factors and is involved in limb induction. Mutations in this gene have been associated with Pfeiffer syndrome, Jackson-Weiss syndrome, Antley-Bixler syndrome, osteoglophonic dysplasia, and autosomal dominant Kallmann syndrome 2. Chromosomal aberrations involving this gene are associated with stem cell myeloproliferative disorder and stem cell leukemia lymphoma syndrome. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]</p>
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Adherens junction, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton

Product images:

Coomassie blue staining of purified FGFR1 protein (Cat# TP320009). The protein was produced from HEK293T cells transfected with FGFR1 cDNA clone (Cat# [RC220009]) using MegaTran 2.0 (Cat# [TT210002]).