

Product datasheet for **TP320030M**

LDLRAD1 (NM_001010978) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human low density lipoprotein receptor class A domain containing 1 (LDLRAD1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220030 protein sequence Red =Cloning site Green =Tags(s)

MNKVFPQGENGYTAAESKAHPGGEAGGGHLCSSRRGACLSASLLLLLATVAALIALVTILGLPSCTPGAQ
ACITLTNRTGFLCHDQRSCIPASGVCDGVRTCTHGEDEDESLCRDVPQSLPHFLVAHCGDPASWIYSDQK
CDGTNNCGDCSELSPVTVCPGPGWWRCPSTFFKYCDCIPRHLCRDHVQHCS DWSDEYACPGP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	21.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001010978
Locus ID:	388633
UniProt ID:	Q5T700



[View online »](#)

RefSeq Size: 2393
Cytogenetics: 1p32.3
RefSeq ORF: 615
Protein Families: Transmembrane

Product images:

Coomassie blue staining of purified LDLRAD1 protein (Cat# [TP320030]). The protein was produced from HEK293T cells transfected with LDLRAD1 cDNA clone (Cat# [RC220030]) using MegaTran 2.0 (Cat# [TT210002]).