

Product datasheet for **TP320163M**

Galectin 8 (LGALS8) (NM_201543) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human lectin, galactoside-binding, soluble, 8 (LGALS8), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220163 protein sequence Red =Cloning site Green =Tags(s)
	<p>MLSLNNLQNIINYNPVIPFVGITPDQLDPGTLIVIRGHVPSDADRFQVDLQNGSSMKPRADVAFHFNPRFK RAGCIVCNTLINEKWGREIITYDTPFKREKSFEIVIMVLKDKFQVAVNGKHTLLYGHHRIGPEKIDTLGIY GKVNIIHSIGFSFSSDLQSTQASSLELTEISRENVPKSGTPQLRPFPAARLNTMPGPGRTVVVKGEVNANA KSFNVDLLAGKSKDIALHLNPRLNIAFVRNSFLQESWGEEERNITSFPFSPGMYFEMIIYCDVREFKVA VNGVHSLEYKHRFKELSSIDTLEINGDIHLLEVRW</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	35.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_963837
Locus ID:	3964




[View online »](#)

UniProt ID:	O00214
RefSeq Size:	6155
Cytogenetics:	1q43
RefSeq ORF:	948
Synonyms:	Gal-8; PCTA-1; PCTA1; Po66-CBP

Summary: This gene encodes a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified LGALS8 protein (Cat# [TP320163]). The protein was produced from HEK293T cells transfected with LGALS8 cDNA clone (Cat# [RC220163]) using MegaTran 2.0 (Cat# [TT210002]).