

## Product datasheet for **TP320225**

### **MAGP1 (MFAP2) (NM\_002403) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human microfibrillar-associated protein 2 (MFAP2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220225 representing NM_002403 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MRAAYLFLFLPAGLLAQGQYDLDPFPDPDHVQYTHYSDQIDNPDIYDYQEVTPRPSEEQFQFQSQQQV QQEVIPAPTPEPGNAELEPTEPGPLDCREEQYPCTRLYSIHRPCKQCLNEVCFYSLRRVYINKEICVRT VCAHEELLRADLCRDKFSKCGVMASGLCQSVAAASCARSCGSC  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	18.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_002394</a>
Locus ID:	4237
UniProt ID:	<a href="#">P55001</a> , <a href="#">A0A024RA94</a>



[View online »](#)

RefSeq Size: 1066

Cytogenetics: 1p36.13


RefSeq ORF: 549

Synonyms: MAGP; MAGP-1; MAGP1

**Summary:** Microfibrillar-associated protein 2 is a major antigen of elastin-associated microfibrils and a candidate for involvement in the etiology of inherited connective tissue diseases. Four transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Sep 2008]

**Protein Families:** Secreted Protein

### Product images:



Coomassie blue staining of purified MFAP2 protein (Cat# TP320225). The protein was produced from HEK293T cells transfected with MFAP2 cDNA clone (Cat# [RC220225]) using MegaTran 2.0 (Cat# [TT210002]).