

Product datasheet for TP320398

Liver Carboxylesterase 1 (CES1) (NM_001266) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human carboxylesterase 1 (monocyte/macrophage serine esterase 1) (CES1), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220398 representing NM_001266 <div> <div>Red</div>=Cloning site <div>Green</div>=Tags(s) </div> <p> MWLRAFILATLSASAAWGHPSPPVVDTVHGKVLGKFVSLEGFAQPVAIFLGIPFAKPPPLGPLRFTPPQP AEPWSFVKNATSYPPMCTQDPKAGQLLSEFTNRKENIPLKLSLCLYLNITPADLTCKNRLPVMVWIH GGGLMVGAASYDGLALAAHENVVVTIQYRLGIWGGFFSTGDEHSRGNWGHLDQVAALRWVQDNIAF GG NPGSVTIFGESAGGESVSVLVLSPLAKNLFHRAISESGVALTSVLVKKGDVKPLAEQIAITAGCKTTTSA VMVHCLRQKTEEELETTLMKFLSLDLQGDPRSQPLLGTVIDGMMLLKTPEELQAERNFHTVPYMGVI NKQEFGWLIPMLMSYPLSEGQLDQKTAMSLWKSYPVLCIAKELIPEATEKYLGGTDDTVKKKDLFLDLI ADVMFVPSVIVARNHRDAGAPTYMYEFQYRPSFSSDMKPKTVIGDHGDELFSVFGAPFLKEGASEEEIR LSKMVMKFWANFARNGNPNGEGLPHWPEYNQKEGYLQIGANTQAAQKLKDKVAFWTNLFACKAVEK PPQ TEHIEL <div> <div>TR</div> <div>TRPLEQKLISEEDLAANDILDYKDDDDKV</div> </div> </p>
Tag:	C-Myc/DDK
Predicted MW:	60.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001257
Locus ID:	1066
UniProt ID:	P23141
RefSeq Size:	2021
Cytogenetics:	16q12.2
RefSeq ORF:	1698
Synonyms:	ACAT; CE-1; CEH; CES2; hCE-1; HMSE; HMSE1; PCE-1; REH; SES1; TGH
Summary:	This gene encodes a member of the carboxylesterase large family. The family members are responsible for the hydrolysis or transesterification of various xenobiotics, such as cocaine and heroin, and endogenous substrates with ester, thioester, or amide bonds. They may participate in fatty acyl and cholesterol ester metabolism, and may play a role in the blood-brain barrier system. This enzyme is the major liver enzyme and functions in liver drug clearance. Mutations of this gene cause carboxylesterase 1 deficiency. Three transcript variants encoding three different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]
Protein Families:	Druggable Genome
Protein Pathways:	Drug metabolism - other enzymes

Product images:



Coomassie blue staining of purified CES1 protein (Cat# TP320398). The protein was produced from HEK293T cells transfected with CES1 cDNA clone (Cat# [RC220398]) using MegaTran 2.0 (Cat# [TT210002]).