

Product datasheet for **TP320401**

FBXL18 (NM_024963) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human F-box and leucine-rich repeat protein 18 (FBXL18), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220401 representing NM_024963 Red =Cloning site Green =Tags(s)

MASSGEDISNDDDDMHPAAAGMADGVHLLGFSDEILLHILSHVPSTDILNVRRTCRLAALCLDKSLIH
TVLLQKDYQASEDKVRQLVKEIGREIQQLSMAGCYWLPGSTVEHVARCRSLVKVNLGCHLTSRLSKML
SALQHLRSLAIDVSPGFDAQLSSECKATLSRVRELKQTLFTPSYGVPCCTSLEKLLLYFEILDRTREG
AILSGQLMVGQSNVPHYQNLRFYARLAPGYINQEVRLYLAVLSDRTPQNLHAFILSVPGSFAESGATK
NLLDSMARNVVDALQLPKSWLNGSSLLQHMKFNPNFYFSRCTLSGGHLIQQVINGGKDLRSLASLNL
SGCVHCLSPDSLLRKAEDDIDSSILETLVASCNLRHLNLSAAHHHSEGLGRHLCQLLARLRHLRSLSL
PVCVADSAPRADRAPAQPAMHAVPRGFGKVRVGVQSCPSPFSGQACQPSSVFWSLKKNLPFLEHLEL
IGSNFSSAMPRNEPAIRNSLPPCSRAQSVGDSEVAIGQLAFLRHLLTAQLPSVLTGSGLVNIGLQCQQQL
RSLSLANLGMGKVVMPALSDMLKHCKRLRDLRLEQPYFSANAQFFQALSQCPSLQRLCLVSRSGTLQP
DAVLAFMARCLQVVMCHLFTGESLATCKSLQSLRSFQAERPALNVVIFPLLHEGLTDVIRDVPLVHLD
EITLFKSRVAEPPNLWW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	78.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_079239](#)

Locus ID: 80028

UniProt ID: [A0A024R857](#), [Q96D16](#), [B4DH63](#)

RefSeq Size: 8253

Cytogenetics: 7p22.1


RefSeq ORF: 2154

Synonyms: Fbl18

Summary: The protein encoded by this gene is a member of a family of proteins that contain an approximately 40-amino acid F-box motif. This motif is important for interaction with SKP1 and for targeting some proteins for degradation. The encoded protein has been shown to control the cellular level of FBXL7, a protein that induces mitotic arrest, by targeting it for polyubiquitylation and proteasomal degradation. Members of the F-box protein family, such as FBXL18, are characterized by an approximately 40-amino acid F-box motif. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains. [provided by RefSeq, Mar 2016]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified FBXL18 protein (Cat# TP320401). The protein was produced from HEK293T cells transfected with FBXL18 cDNA clone (Cat# [RC220401]) using MegaTran 2.0 (Cat# [TT210002]).