

Product datasheet for **TP320433L**

CORO6 (NM_032854) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Homo sapiens coronin 6 (CORO6), 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC220433 representing NM_032854

Red=Cloning site **Green**=Tags(s)

MSRRVVRQSKFRHVFGQAAKADQAYEDIRVSKVTWDSSFCVNPFLAIIVEAGGGGAFIVLPLAKTGRV
DKNYPLVTGHTAPVLDIDWCPHNDNVIASASDDTTIMVWQIPDYTPMRNITEPIITLEGHSKRVGILSWH
PTARNVLLSAGGDNVWVGTGEVLLSLDDMHPDVIHVCWNSNGSLLATTCKDKTLRIIDPRKGQVV
AERFAAHEGMRPMAVFTRQGHIFFTRMSQRELGLWDPNNFEEPVALQEMDTSNGVLLPFYDPDSS
I

VYLCGKGDSSIRYFEITDEPPFVHYLNTFSSKEPQRGMGFMPKRGLDVSKCEIARFYKLHERKCEPIIMT
VPRKSDLFQDDLYPDTGPEPALEADEWLSGQDAEVLISLRDGYVPPKHRELVRTKRNILDVRPPSGPR
RSQSASDAPLSQQHTLETLLEEIKALRERVQAQEQRITALENMLCELVDGTD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 52.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: NP_116243

Locus ID: 84940


UniProt ID: Q6QEF8

RefSeq Size: 2450

Cytogenetics: 17q11.2

RefSeq ORF: 1416

Product images:



Coomassie blue staining of purified CORO6 protein (Cat# [TP320433]). The protein was produced from HEK293T cells transfected with CORO6 cDNA clone (Cat# [RC220433]) using MegaTran 2.0 (Cat# [TT210002]).